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*In-Office dental whitening gel*

**For professional use only**

*Carefully read all information in this document before using the product. Keep it for reference at least until the product has been fully consumed and/or until there is no further interaction between this product and the last patient for whom it will be used.*

**Product Description**

Whiteness HP AutoMixx 6% is a hydrogen peroxide-based whitener for whitening vital teeth, for in-clinic use only. It is presented in a double-barrel syringe with two compartments, one containing hydrogen peroxide and the other containing a thickener. The phases are mixed automatically – by means of a mixing tip – to form the whitening gel, which will be ready to be applied to the teeth. The gel has a blue color for better identification of application areas and to facilitate removal.

The product’s mechanism of action is based on the chemical reaction between hydrogen peroxide and the pigments within the tooth structure (enamel and dentin), and this reaction does not depend on the application of light sources (e.g.: LED, laser, etc.). Thus, the use of light sources is optional and the professional’s responsibility.

**Presentation form**

1 double-barrel syringe with 4 g of whitening gel, 4 mixing tips and 4 brush tips.

**Basic composition**

After mixing the phases: hydrogen peroxide, thickeners, neutralizer, calcium compound, glycol, dye, inorganic filler, purified water and aroma.

**Indications for use**

The product is suitable for the following cases:

- Naturally dark or yellowed vital teeth;
- Vital teeth darkened by age;
- Vital teeth that are darkened or stained by extrinsic substances, such as pigments from coffee, tobacco, tea, etc.;

- Vital teeth with mild to moderate staining due to fluorosis;

- Vital teeth stained by tetracycline, in cases where the staining is mild;

- In more challenging cases (severe darkening), it can be used in combination with the home whitening technique.

**Precautions**

- The product is for clinic use only and should only be handled by dentists;
- The patient must be in good oral and general health before starting treatment;

- When handling the product, both the professional and the assistant must wear protective gloves and goggles. The patient must also wear protective eyewear and other necessary protections to avoid accidental product contact with the skin and clothing;

- Before opening the package, protect your hands and eyes and make sure there is no leaked liquid inside. In case of leakage, request replacement of the material;
- Handle the product away from the patient, over a sink or safe space without risk of contaminating other people;
- After using the product, wash hands and materials contaminated with peroxide thoroughly. Think about prevention and safety when handling the product;

- Patients should be advised to avoid eating acidic and/or strongly colored foods for at least 24 hours after whitening to prevent sensitivity and possible jeopardizing the whitening.

**Contraindications**

- Acid etching of tooth enamel prior to whitening is neither necessary nor recommended;

- The product is not recommended for people who have recently undergone gum surgery or who have periodontal disease, sensitive teeth, exposed roots, cavities, ulcers or general poor health. In such cases, a solution to the problem should be provided before starting treatment;

- The product is not indicated for whitening teeth with imperfect amelogenesis and dentinogenesis, severe fluorosis, intense tetracycline staining or other enamel and dentin anomalies that jeopardize tooth vitality;

- The product is not suited for use in patients under anesthesia. It is important to monitor the sensitivity felt by the patient during whitening. In extreme cases, when there is some imperfection or flaw in the dental structure that makes it possible for peroxide to reach the dental pulp in greater quantity, pulp death can occur;

- The product is not suitable for whitening in pregnant women, lactating women and patients under 12 years of age;

- The product is not recommended for people with allergic reactions to the product’s components;
- The product is not suitable for use with high-power lasers and others that deliver high temperatures to the dental surface. The pulp temperature should not exceed 42°C (108°F), to prevent irreversible damage to its vitality.

**Adverse effects**

Dental hypersensitivity can occur during or after whitening. As this is often a difficult condition to predict, using Desensibilize KF 0.2% (FGM) – for home use – or Desensibilize KF 2% (FGM) – clinic use – is recommended as prophylactic desensitization before starting whitening.

Tooth whitening can increase the visibility of pre-existing enamel hypoplasia (opaque white patches). In most cases, no additional care is needed because the stains disappear as soon as the whitening is finished and the tooth resumes its usual hydration. For more resistant cases, treatment with fluoride in the highlighted region can help it disappear.

**Directions for use**

Preparing for whitening:

- Suit up (dentist and assistant) with your safety equipment (gloves, goggles, cap, mask and apron) to handle the product. Also equip your patient to prevent accidental contact with the gel, using a cap, goggles and a suitable lip retractor (Arcflex – FGM).

- Remove the whitener from the packaging observing the safety recommendations.

- Away from the patient, over a sink or counter, attach the mixing tip (make sure the mixing tip fits the syringe perfectly).

- Press the plunger until the product comes out of the tip in gel form. Small amounts of unmixed product may be found at the beginning of mixing and should be discarded.

**Note:** If the whitener is stored under refrigeration, remove it at least 30 minutes before starting the treatment so that its temperature matches that of the environment.

**Whitening of vital teeth**

Perform a good assessment of the patient’s oral cavity (read precautions and contraindications).

- Prophylaxis is recommended.

**Step-by-step instructions for vital teeth**

- Select and record the color of the patient’s teeth using a color chart and/or photograph before starting whitening.

- After recording the color, insert the Arcflex lip retractor (FGM) to facilitate access and manipulation of the oral cavity. Then apply Desensibilize KF 2% (FGM) to the teeth to be whitened for 10 minutes before starting whitening.

**Note:** Due to the low concentration of AutoMixx 6% whitening gel, prior application of a desensitizer is optional.

- Remove the desensitizing gel with water spray and dry the gum tissue. Establish the relative isolation with Top dam (light-curing gingival barrier - FGM), covering the marginal gum and papillae with a layer 3 to 5 mm wide and no more than 1 mm thick. The barrier should cover approximately 0.5 to 1 mm of the tooth surface and gum to ensure a good seal. Use a clinical mirror, looking from incisal to cervical, and check for uncovered gingival tissue. If there is, correct it. After proper

application, light cure the Top dam composite (FGM) for 20 to 30 seconds for each group of 3 teeth. The gingival protector that forms is rigid and insoluble, preventing possible irritation by the whitening agent.

**Note:** Due to the low concentration of AutoMixx 6% whitening gel, use of a light-curing gingival barrier during the procedure is optional.

- Press the plunger gently until the two phases (peroxide and thickener) slowly mix in the mixing tip. A small amount of gel should be discarded in a suitable container, ensuring that the applied product is properly blended.

- Apply the gel directly to the teeth and surfaces to be whitened (usually from the 2nd premolar to the other on the vestibular/buccal surface). Apply the finest layer of gel possible, using the tip itself to spread the product.

- Once the gel is applied to the teeth, it should be left in contact for up to 50 minutes. During application, air bubbles that are formed by the gel action should be removed to renew contact of the gel with the teeth. There is no need to apply a new layer of gel during the session.

- At the end of the treatment, suction out the gel and wash the teeth. When used, remove the gingival protector by detaching it with an exploratory probe.

- Polish the teeth with Diamond Excel Polishing Paste (FGM) and Diamond (FGM) or Diamond Flex (FGM) felt discs.

9. A return visit should be schedule after a 2- to 3-day break, allowing for effective evaluation of the color of teeth. If necessary, additional whitening sessions can be performed following the steps described above. Generally, 2 or 3 sessions are enough.

**Additional information**

**General:**

- Whiteness HP AutoMixx 6% can be used in combination with supervised at-home whitening. According to this technique, 2 or 3 sessions of Whiteness HP AutoMixx 6% are performed in the dental clinic, with a 2- to 3-day break between each session. In parallel, supervised home bleaching is performed for 3 to 4 weeks;

- The results depend on the origin and age of the stains. In some cases, more than 3 in-clinic sessions may be required to achieve satisfactory results;

- If the patient develops hypersensitivity that cannot be controlled, the professional will need to evaluate whether there is any local factor that may be causing the discomfort. If there is no apparent cause, consider reducing the application times in the clinic. If this hypersensitivity cannot be controlled, stop the treatment;

- As with other whitening techniques, we recommend recording the teeth color before whitening, informing the patient about the real possibilities in his/her case (limitations in cases of tetracycline staining, gray coloration, etc.), warning about possible hypersensitivity and the need to change restorations. There are cases in which, due to particular characteristics of the teeth (enamel characteristics, type of pigmentation, etc.), the desired level of whitening may not be possible.

**About the gel:**

- Right after mixing, the product has a bluish color that does not change during whitening;

- The product was designed so that each double-bore syringe containing 4 g of product is sufficient for 3 or 4 complete sessions (each session corresponding to application on the upper and lower arches, from the

2nd premolar on one side to the other). However, this can change depending on the amount of gel applied to the surface of teeth. For best efficiency, apply the finest possible layer;

- When gel remains in the syringe, it can be stored (see storage requirements) for later use until it expires (see packaging).

**Note:** The decision to use any light source on the product is the professional’s sole and total responsibility.

**Storage conditions**

- Keep the syringe with its original, properly adapted cap;

- Store the product at temperatures between 5°C and 25°C (41°F and 77°F) for best conservation;

- Do not freeze the product;

- Protect it from direct sunlight.

**Handling conditions**

Mixing tips are disposable and should be disposed of properly after use.

For product disposal, follow the laws in your country.

**Warnings**

- Do not use the product if its expiration date has passed;
- Keep out of the reach of children;
- Irritating to eyes and skin.

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